

Samplename : Adv\_2  
Studyowner : walkup  
DayQ : PROTON  
NightQ

---

Date: 2025-11-03 at 13:41:03, Adv\_2 Submitted  
Date: 2025-11-03 at 14:07:11, Adv\_2\_20251103\_01 Created  
Date: 2025-11-03 at 14:07:12, Adv\_2\_20251103\_01 dayQ Started  
Date: 2025-11-03 at 14:07:13, Adv\_2\_20251103\_01: PROTON\_001 Setup  
Date: 2025-11-03 at 14:07:15, preScan: Pre-acquisition delay Start  
Date: 2025-11-03 at 14:08:48, preScan: Pre-acquisition delay done  
Date: 2025-11-03 at 14:08:49, preScan: Lock sample Start  
Date: 2025-11-03 at 14:08:57, preScan: Lock sample done  
Date: 2025-11-03 at 14:08:58, preScan: Probe tuning Start  
Date: 2025-11-03 at 14:09:03, preScan: Probe tuning Start  
Date: 2025-11-03 at 14:09:09, preScan: Probe tuning Start  
Date: 2025-11-03 at 14:09:24, preScan: Probe tuning Start  
Date: 2025-11-03 at 14:09:37, preScan: Probe tuning done  
Date: 2025-11-03 at 14:09:39, preScan: Gradient shimming Start  
Date: 2025-11-03 at 14:10:19, preScan: Gradient shimming done  
Date: 2025-11-03 at 14:10:19, preScan: Auto gain Start  
Date: 2025-11-03 at 14:10:26, preScan: Auto gain Start  
Date: 2025-11-03 at 14:10:33, preScan: Auto gain done  
Date: 2025-11-03 at 14:10:33, Adv\_2\_20251103\_01: PROTON\_001 Acquisition started  
Date: 2025-11-03 at 14:13:30, Adv\_2\_20251103\_01: PROTON\_001 (Error) Auto lock fa  
ilure in finding resonance of sample